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
www.baumann-dental.de

2017 - 08 

3D AESTHETIC
FACEBOW BRIEF
INSTRUCTIONS

ENGLISH

The basic functionality of **the 3D aesthetic facebow** does not differ from our classical Arto facebow (#72100). We therefore refer expressly to its instructions for use, which are included in the scope of delivery of the 3D aesthetic facebow. Should the instructions be missing, or should you not understand them completely, please contact our Customer Service before using this product:
Phone +49(0)7236-933690

The 3D aesthetic facebow complements the basic functionality of the standard Arto facebow with the option of correcting the reference plane through three axes. Longitudinal axis = the bipupillary line is applied as horizontal reference from the front. The independently height-adjustable earpieces allow adjustment of the longitudinal axis on the patient as the human auditory canals are not absolutely transversal.

Transverse axis = the "reference plane" can be adjusted via the height-adjustable glabella. Its scale describes the established reference planes as "mean value": Camper's plane (Campersche Ebene, CE), patient's plane (Patienten-Ebene, PE) or the Frankfurt horizontal (Frankfurter Horizontale, FH). Infinite adjustment allows for precision alignment on the patient. The plane indicator facilitates finding the respective reference points on the cranium.


Vertical axis = the sagittal plane passes through the centre of the cranium and corresponds with the palatal suture (from a caudal view). Precise alignment of the sagittal plane is enabled by the infinite lateral movement of the glabella.

Adjustment of the sagittal plane

Adjustment of the sagittal plane on the patient by aligning the facebow to the centre of the face. The centre of the facebow acts as reference (screw channel of the joint support).

Registration of the transversal


Setting the bipupillary plane as horizontal reference from the front by adjusting the height-adjustable earpieces.

Attention: to minimise measuring errors, set all settings to the scale value "0" prior to every application as a matter of principle. This refers to the height-adjustable earpieces and the lateral adjustment of the glabella. Then set the desired reference plane.

Recommendation: registration should preferably be performed on patients sitting in an upright position. Ideally, the assistant stands behind the patient, fixates the head and supports the clinician in applying the bite fork and bow.

Comparison with
"natural head posture"


The attachable spirit level allows a comparison of the adjusted longitudinal axis on the cranium (bipillary plane) with the patient's "natural head posture". The spirit level indicates two axial directions, these being the longitudinal and the transverse axes.

CE - Camper's plane

The illustration shows the alignment of the facebow to the Camper's plane (CE). The height-adjustable glabella is set to "CE", whereby the plane indicator is set to the subnasal point (spina nasalis anterior).

CE - Camper's plane


The height-adjustable glabella is aligned along the Camper's plane as reference plane. The plane indicator points to the spina nasalis anterior.

Cast mounting in the articulator: direct coupling of the 3D aesthetic facebow to the articulator is not intended. For information on the transfer of the registration to the articulator, please refer to the overall instructions for our Arto facebow (#72100).

Recommendation for setting the articulator: if no information on the patient's condylar path inclination is available, the following settings are recommended for the articulator, for registration at the reference plane:

- Camper's plane (CE): 30° condylar path inclination
- Patient's plane (PE): 35° condylar path inclination
- Frankfurt horizontal (FH): 45° condylar path inclination

PE - patient's plane


The illustration shows the alignment of the facebow to the patient's plane. The height-adjustable glabella is set to "PE". The plane indicator is set to the nasal wing.

The 3D aesthetic facebow was developed in cooperation with Dr. med. habil. Volker Ulrici and DMT Walter Böthel from Leipzig. Always follow the current legal requirements for the hygienic processing of medical devices! This relates to the following components of the facebow: earpieces, Nasion adapter, bite fork.

Recommendation: the Robert Koch Institute (www.rki.de) provides information on the processing of medical devices on its website under the section on "Infection protection".

FH - Frankfurt horizontal

The illustration shows the alignment of the facebow to the Frankfurt horizontal. The height-adjustable glabella is set to "FH".

FH - Frankfurt horizontal


The illustration shows the alignment of the facebow to the Frankfurt horizontal. The plane indicator is set to the infra-orbital point.

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Baumann Dental GmbH